

Process Diagnostics for Pulp and Paper

Moisture Balance



Moisture Balance

Description

The Moisture Balance100 (MB100) is a rapid, repeatable weight loss moisture balance specifically tuned to optimize MasterScreen sample testing. Obtaining the consistency of sample to be tested is required to generate repeatable contaminant information. The MB100 generates rapid, oven dried weights for filter paper and pulp samples, in consort with a speed dryer. If gramametric results are desired, the MB100 is used to generate the oven dried weight results on the rejects and or/accepts generated by the MasterScreen.

Benefits

The MB100 reduces tester contact time and increases the accuracy and repeatability of generating moisture free weight data.

Features

- Customized programmable settings
- · Easy to use and install
- Easy to read display screens
- · Provides rapid, oven dried weights for testing

Space and Services Required

- 120 V, 50 Hz
- Dimensions: $13 \times 9 \times 7$ in (W x L x H)
- Weight: 15 lbs

Pulmac Systems International Inc. 35 Dorset Lane, Suite 201 Williston, Vermont 05495 USA

> Phone: + 802.223.2663 Fax: + 802.223.2664 Email: info@pulmac.com URL: www.pulmac.com